



Product Introduction

BeaverBeadsTM Mag NH₂magnetic beads have following characteristics: superparamagnetism, rapid magnetic responsiveness, rich amino functional groups, monodispersity and sub-micron scale size. When it is in special chemical reagent such as glutaraldehyde, polypeptide, protein, oligonucleotide and other bioligands can be covalent to the surface of microspheres. BeaverBeadsTM NH₂magnetic beads are important carriers of the medical and molecular biology research tools

BeaverBeads[™] Magrose NH₂ magnetic beads adopt advanced polymer technology to combine super paramagnetic materials and polymer materials to form a new functional magnetic microsphere. Compared with the traditional magnetic beads, Magrose NH₂ magnetic beads have faster magnetic responsiveness, while keeping the microspheres good dispersibility, low non-specific adsorption and richer binding spots, can easily and efficiently reach high load combination with a variety of biological ligands (proteins, peptides, oligomeric nucleotide, drug molecules, etc.), can be used as good material for covering, adsorption, chemical modification etc. subsequent processing.

Product Information

Product name	Mag NH ₂ -500	Mag NH ₂	Magrose NH ₂		
Mean particle size	500nm single dispersion*	2 µm single dispersion	30~150 μm		
Surface carboxyl/content	NH ₂ (~40 μM/g)	NH ₂ (~200 μM/g)	NH ₂ (~50 μM/mL gel)		
Magnetic core	Fe ₃ O ₄	Fe ₃ O ₄	Fe ₃ O ₄		
Shell	Silica	Silica	Magrose		
Magnetic	superparamagnetism	superparamagnetism	superparamagnetism		
Saturation magnetization	50.11 emu/g	45.23 emu/g	41.09 emu/g		
Specific Surface Area	20.01 m ² /g	7.16 m ² /g	/		
Preservation solution		2 years, store at 2 °C~8°C (long term preservation)			
* water average particle size, Malvern Nano determination					

Product Advantages

- 1. Rich binding spots, more specific binding with ligands.
- 2. Superparamagnetism and high magnetic responsiveness, saving operation time.
- 3. Good dispersion and resuspension ability, improving the operation process.
- 4. Good physical and chemical stability, ensuring repetitive effects.

Coupling method of magnetic beads and biomolecules (reference, taking protein A as an example)

1. Pre-treat the magnetic beads

After blending magnetic beads, take 100µL Mag/Magrose magnetic beads into 1 mL EP tube, magnetic separation, then remove the supernatant, then washing three times with 200µL PBS solution (50 mM PBS, pH 7.4), magnetic separation and discard the supernatant;

2. Activate Glutaraldehyde

Add fresh 100 μ L glutaraldehyde solution (15%) to the EP tube, vortex mix magnetic beads to fully suspend, cover with tin foil paper for 1 h at 25 °C to react (avoid lights to avoid glutaraldehyde aggregation). During the reaction, keep the magnetic beads suspended (vertical mixing apparatus can help);

3. Wash after the activation

古非命科察阳弱

After the magnetic beads were activated, magnetic separation and discard the supernatant, and then wash three times with 50 mM PBS and pH 7.4 buffer.

www.beaverbio.com

400-600-3979

4. Couple with magnetic beads

Add 50µL ~200µL biological ligands to the above EP tube containing magnetic beads (appropriate volume and concentration needs to be optimized, keep the solution pH 8.0, add 0.05% Tween - 20 when necessary in order to improve the magnetic beads dispersion), gently blending; cover with tin foil paper to couple 3 h (dark environment) at 25 °C, or 25 °C coupling for 1 h and after that place at 4 °C coupling over the night, keep magnetic beads suspended when coupling (vertical mixing apparatus can help);

Service@beaverbio.com

Note: the glutaraldehyde solution at 280 nm has absorption peak, so the coupling protein content with the magnetic beads cannot be calculated by OD280 value difference before and after the reaction. The content can be indirect measured by BCA kit or protein electrophoresis method.

5. Close after coupling

Place EP tubes in the magnetic separator, remove the supernatant, add 200 ~ 1000µL BSA/PBS solution (pH 7.2, containing 5% BSA) to resuspend the magnetic beads (use ultrasonic when necessary), react for 1 h at 25 °C to close magnetic bead surface nonspecific adsorption spots, keep the magnetic beads suspended (vertical mixing apparatus can help);

6. Store

Place EP tubes in the magnetic separator, remove the supernatant with 200µL PBS solution (pH 7.2) or wash three times by storage solution, resuspend in the suspension (the storage solution amount are determined by the concentration of the coupling ligand magnetic beads), and finally store at 4 °C. When necessary, add 0.02% (w/v) sodium azide (NaN3) as the preservative to inhibit bacterial growth.

Note

1. This product should not be frozen, dried or centrifuged. Freezing, drying and centrifugation will cause the beads to agglomerate, not easy to resuspend and disperse, and affect the chemical activity of beads surface functional group.

Before using this product, be sure to fully oscillate or ultrasonic to keep the beads in a uniform suspension.
In the process of use, to remove the ethanol in the storage solution, use pure water or buffer to wash magnetic beads 2~3 times.

4. This product must be used with magnetic separation equipment.

5. In order to ensure optimal experimental results, please select appropriate ligands for covalent coupling reactions.

6. This product is for research use only.

Product list

NO.	Product name	Specification	Size	Concentration
70201-5	BeaverBeads TM Mag $NH_2 - 500$	5mL	500 nm	10 mg/mL
70201-50	BeaverBeads [™] Mag NH ₂ – 500	50mL	500 nm	10 mg/mL
70202-5	BeaverBeads [™] Mag NH ₂	5mL	2 µm	10 mg/mL
70202-50	BeaverBeads [™] Mag NH ₂	50mL	2 µm	10 mg/mL
70203-5	BeaverBeads TM Magrose NH_2	5 mL	30~150 µm	20%(v/v)
70203-50	BeaverBeads [™] Magrose NH ₂	5mL	30~150 µm	20%(v/v)

Limited use of trademark license

Beaver Biomedical Engineering Co., Ltd. disclaims all intellectual property rights owned or co-owned by the developer and all the contents and services developed or cooperated with other units and is protected by intellectual property laws such as trademark, patent and copyright. The rights of the purchaser of this product are limited to the use of the purchaser of the product for internal research and the right is not transferable and is not available for any commercial application and the purchaser is not entitled to re-sell the product or any part thereof The Such as for commercial use (including but not limited to sales of the agent), must be approved by Beaver Biomedical Engineering Co., Ltd., and when used to indicate the source and intellectual property rights, copyright, etc. Beaver Biomedical Engineering Limited Company all the mark. For more information, please contact Beaver@beaverbio.com or Beaver Biomedical Engineering Co., Ltd. Address: Biobay A6-101, 218 Xinghu Street, Suzhou Industrial Park. This product by Beaver Biomedical Engineering Co., Ltd.





www.beaverbio.com

Service@beaverbio.com

400-600-3979

Copyright Notice:

© 2013 Beaver Biomedical Engineering Co., Ltd. All rights reserved. Any of the contents of this User's Guide, regardless of trademark, design, text, images and any other information without special instructions, are copyrighted by Beaver Biomedical Engineering Co., Ltd. . The Company reserves the right to take legal action and to hold its duties in violation of the relevant laws and regulations of the State and does not respect the statement, without the consent and unauthorized use of the contents of this User Manual.

For support, please visit: www.beaverbio.com/support or email: Service@beaverbio.com