

## BeaverBeads™ Streptavidin

### Product Introduction

(SA-Biotin) system has extremely high binding affinity ( $K_d=10^{-15}$ ). It has widely applications in the biological field. BeaverBeads™ Streptavidin covalently connects the SA to the solid carrier surface by using the patented protein coupling technique of Beaver, which can efficiently bind the ligands such as biotinylated antibodies, nucleic acids, proteins etc. The product adopts superparamagnetic microspheres with uniform size and regular morphology, which facilitates the rapid and convenient capture of target molecules and the realization of magnetic separation. This product can be equipped with automation equipment for high flux operation.

### Product Information

Product Information	BeaverBeads™ Streptavidin (300nm)
<b>Biotinylated IgG Binding Capacity</b>	> 2 µg/mg beads
<b>Bead Concentration</b>	10 mg/mL
<b>Bead Surface</b>	hydrophilic group
<b>Preservative Solution</b>	1xPBS(10mM PBS) pH 7.4, including 0.1% (v/v) Tween -20, 0.1% (w/v) NaN <sub>3</sub>
<b>Storage Condition</b>	2~8°C
<b>Shelf Life</b>	12 months

### Product Application Area

Suitable for IVD industry.

### Chemiluminescence Immunodiagnostic Process

(this process takes a two-step double antibody sandwich CLIA as an example)

#### 1. Preparation

- Reagents: capture antibody, analyte standard, enzyme-labeled antibody, substrate solution, etc., adjust to room temperature before using.
- Washing Buffer: Adjust it to room temperature before using.
- Chemiluminescent 96-well flat bottom microtiter.
- 96-well flat bottom magnetic separator: Beaver magnetic separator, Cat. No.60302.
- Vortex generator.
- Vacuum suction pump.
- Chemiluminometer
- Pipette.

#### 2. Magnetic Particle Chemiluminescence Immunodiagnostic Process

- 1. Make sure that the beads have been adjusted to the proper concentration(0.8 mg/mL). Place the beads on the vortex generator for 20s and resuspend the beads with oscillation. Pipette 50 µL magnetic beads of the appropriate concentration into a 96-well plate, magnetic separation, discard the supernatant with a pipette, remove 96-well plate from the magnetic separator.
- 2. Add 100 µL biotinylation capture antibody, fully resuspend the magnetic beads, incubate at 37°C for 15 min, magnetic separation, discard the supernatant with a pipette, and remove 96-well plate from the magnetic separator.
- 3. Add 200 µL Washing buffer in each well, fully resuspend the magnetic beads, magnetic separation, discard the supernatant with a pipette, and remove the 96-well plate from the magnetic separator, and

repeat this step 2 times, totally wash 3 times.

- 4. Add 50 µL tested standard sample or tested sample in each well, fully resuspend the magnetic beads, incubate at 37°C for 15 min, magnetic separation, discard the supernatant with a pipette and remove 96-well plate from the magnetic separator.
- 5. Add 200 µL Washing buffer in each well, fully resuspend the magnetic beads, magnetic separation, discard the supernatant with a pipette, and remove the 96-well plate from the magnetic separator, and repeat this step 2 times, totally wash 3 times.
- 6. Add 100 µL enzyme-labeled antibody, fully resuspend the magnetic beads, incubate at 37°C incubator for 15 min, magnetic separation, discard the supernatant with a pipette, remove 96-well plate from the magnetic separator.
- 7. Add 200 µL Washing Buffer in each well, fully resuspend the magnetic beads, magnetic separation, discard the supernatant with a pipette, remove the 96-well plate from the magnetic separator, and repeat this step 2 times, totally wash 3 times.
- 8. Add 150 µL substrate solution in each well, fully resuspend the magnetic beads, and incubate for 5 min in the dark.
- 9. Place the 96-well plate into the chemiluminometer and perform the appropriate data processing.

### Note

1. Avoid drying, freezing and high-speed centrifugation of beads
2. In order to reduce the loss of magnetic beads, the time of magnetic separation should be no less than 1 min.
3. Recommend to use a low-adsorption pipette tip and reaction tube to avoid loss due to adherent beads and solutions.
4. Binding capacity of magnetic beads and biotinylated molecules are related to the size of the molecular weight. Users need to determine the loading of magnetic beads on specific biotinylated molecules according to experiments.
5. This operation process is suitable for manual chemiluminescence immunoassay. If using automatic or semi-automated chemiluminometers operation, please refer to the instrument manual.
6. This product is for research and use only.

### Product List

Cat. No.	Product Name	Specification
22308-1	BeaverBeads™ Streptavidin	300 nm, 1 mL, 10 mg/mL
22308-10	BeaverBeads™ Streptavidin	300 nm, 10 mL, 10 mg/mL
22308-100	BeaverBeads™ Streptavidin	300 nm, 100 mL, 10 mg/mL
60201	Magnetic Separator Stand 2/15	1/Pk., suitable for 1.5 mL, 2 mL EP tubes and 15 mL centrifuge tubes
60203	Magnetic Separator Stand 50	1/Pk., suitable for 50 mL centrifuge tubes
60304	Magnetic Separator Stand 96 III	1/Pk., suitable for 96 hole Deep-well Multiwell Plate

#### Limited use of trademark license

Beaver Biomedical Engineering Co., Ltd. disclaims all intellectual property rights owned or co-owned by the developer and all the contents and services developed or cooperated with other units and is protected by intellectual property laws such as trademark, patent and copyright. The rights of the purchaser of this product are limited to the use of the purchased quantity of the product for internal research and the right is not transferable and is not available for any commercial application and the purchaser is not entitled to re-sell the product or any part thereof. The Such as for commercial use (including but not limited to sales of the agent), must be approved by Beaver Biomedical Engineering Co., Ltd., and when used to indicate the source and intellectual property rights, copyright, etc. Beaver Biomedical Engineering Limited Company all the mark. For more information, please contact [Beaver@beaverbio.com](mailto:Beaver@beaverbio.com) or Beaver Biomedical Engineering Co., Ltd. Address: Biobay A6-101, 218 Xinghu Street, Suzhou Industrial Park.  
This product is produced by Beaver Biomedical Engineering Co., Ltd.

#### Copyright Notice:

© 2013 Beaver Biomedical Engineering Co., Ltd. All rights reserved. Any of the contents of this User's Guide, regardless of trademark, design, text, images and any other information without special instructions, are copyrighted by Beaver Biomedical Engineering Co., Ltd. The Company reserves the right to take legal action and to hold its duties in violation of the relevant laws and regulations of the State and does not respect the statement, without the consent and unauthorized use of the contents of this User Manual.

For support, please visit: [www.beaverbio.com/support](http://www.beaverbio.com/support) or email: [Service@beaverbio.com](mailto:Service@beaverbio.com)